



YOUR VIRTUAL TREATMENT STARTS HERE




CAPTURE YOUR SMILE




SHARE IT WITH THE DOCTOR



RECEIVE FEEDBACK



SEE YOUR PROGRESS



DOWNLOAD THE APP



WHEN YOU START TREATMENT WITH OUR OFFICE, YOU'LL RECEIVE AN EMAIL FROM DENTAL MONITORING

1. Open the email and click the download link.
2. Download the app and login with your email and password provided from the email.
3. Click 'Enable Notifications' when prompted - this is how you'll be alerted to take your next scan.
4. Practice! View the tutorial to familiarize yourself with taking your first scan. FAQs can be found in the help section in the top right hand corner.

ALIGNER DELIVERY



WHEN YOU RETURN TO OUR OFFICE TO RECEIVE YOUR ALIGNERS, MAKE SURE TO BRING YOUR SMARTPHONE

1. Open the DM app.
2. Click 'Take Scan.'
3. Insert your cheek retractors and follow the prompts on the screen
4. Scan and show off your pearly whites! They shouldn't be covered by your lips, tongue or the retractor. Avoid backlight, and find neutral, plain backgrounds.

SCAN!



WHEN A NEW SCAN IS REQUIRED, YOU WILL RECEIVE A NOTIFICATION ON YOUR SMARTPHONE

1. Open the DM app.
2. Click 'Take Scan.'
3. You will receive a notification within 12 hours if you can move on to your next aligner, or remain on your current set for a few more days.
4. If you have any issues with the app, please email support@dental-monitoring.com