

Carriere Motion Appliance



The Carriere Motion Appliance is shown with Class II elastic attached to hook on the cuspid pad.

- Used at the beginning of treatment when the patient is most motivated to achieve results. Good results in this phase can significantly reduce the amount of time spent in braces.
- Wear your retainer and elastics 22+ hours per day to achieve tooth movement.
- The Carriere Motion Appliance corrects the bite by rotating and uprighting the upper first molars while moving the canine and premolars backward.
- Independently moves each posterior segment, from canine or premolar to molar, as a unit.
- Please see instructions sheets on retainer care and recommended foods to avoid for additional information